



Next Generation Kits For DNA Extraction

Pure DNA with only one centrifugation step and no Bind-Wash-Elute requirements

Ideal for use in teaching Labs and basic research

Back to the Lab - 50% OFF

Description	Code	Price	Sale
EchoLUTION Blood DNA Micro Kit (50)	BIOE010-001-050	\$270	\$135
EchoLUTION Cell Culture DNA Kit (50)	BIOE010-006-050	\$310	\$155
EchoLUTION Plant DNA Kit (50)	BIOE010-003-050	\$340	\$170
EchoLUTION Tissue DNA Micro Kit (50)	BIOE010-002-050	\$350	\$175
EchoLUTION Viral RNA/DNA Swab Kit Plus (50)	BIOE012-002-050	\$420	\$210



Based on the principle of reverse chromatography, BioEcho extraction kits use a specialised matrix allowing DNA fragments to flow through the column with unwanted sample-derived contaminants being bound.

[more details](#)



BioEcho is an environmentally friendly science providing no harsh chemicals, recyclable containers and providing 70% less plastic than a standard silica kit.

VALID UNTIL:
30 Apr 2023
Ref. Code:
BIOE-DNA50