



## Dual Chamber -30°C Biomedical Freezer

2 completely independent chambers with their own temperature controls

Description	Code	Price	Sale
PHCbi Upright Dual Chamber -30°C Biomedical Freezer	SANYMDF-MU539D-PE	\$9,841	<b>\$7,995</b>

- 479L capacity (approx. 240L per chamber)
- Supplied complete with 6 shelves & 12 bins: store in bins or on adjustable height shelves
- Temperature range: selectable from -20°C to -30°C in each chamber
- Push button manual defrost for each chamber
- Separate door for each chamber
- Direct cooling
- Drainage hose for easy removal of defrosted water
- Alarms for hi & low temperature and power failure
- Lockable doors for security
- Quiet: only 42dB(A)
- Filter-less design – no messy filter(s) to remove/clean/reinstall
- Access Ports for remote probe entry
- 3-year warranty



VALID UNTIL:  
30 Apr 2023  
Ref. Code:  
SANY-Dual

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

