

TwinGuard 360L ULT Freezer

Their most secure ultra-low temperature freezer
for the storage of your valuable samples

Code	Description	Price	Sale
SANYMDF-DU302VX-PE	360L -86°CULT Upright Freezer	\$21,911	\$14,895

Only While Stocks Last!




TwinGuard ULT freezers have two independent refrigeration systems to provide the highest levels of sample security.

PHCbi's Dual Cooling System's backup system ensures that the integrity of irreplaceable and potentially life-saving samples are not put at risk.

Two individual compressors provide a reliable -86°C Ultra-Low temperature environment. Should one cooling circuit unexpected fail the other will maintain a temperature of at least -70°C.

- Temperature range: adjustable from -50°C to -86°C
- Filter-less design – no messy filter to clean!
- VIP Plus Vacuum Insulation Panels
- 5-Year Warranty



VALID UNTIL:
31 May 2023
Ref. Code:
SANY-ULTTwin

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.



Compact Laboratory -30°C Chest Freezer 138L

Only While Stocks Last!

Description	Code	Price	Sale
PHCbi Chest Freezer -30°C 138L	SANYMDF-137-PE	\$4,059	\$2,750

Includes 2 baskets for sample storage

A storage environment with excellent safety features, easy operability, and a host of other features, these freezers offer unsurpassed reliability and functionality.



- Push button manual defrost
- Quiet in operation: 35dB(A)
- Touchpad control panel with digital temperature display
- Lid lock for security
- Drain port for removal of water from bottom of chamber
- Alarms for high temperature, low temperature, and power failure
- Dry remote alarm contacts for connection to BMS
- 2 access ports as standard
- Castors for mobility and locating feet when in working position
- 5-year warranty

VALID UNTIL:
31 May 2023
Ref. Code:
SANY--30BIOM

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in AUD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.



Large Capacity Pharmaceutical/Vaccine Refrigerators

FREE remote alarm/power failure alarm
(valued at \$382)

These units deliver stable, reliable temperature control for clinical, pharmaceutical and biomedical storage.

Code	Description	Price	Sale
SANYMPR-1412-PE	1,364L with 4 adjustable shelving levels (8 shelves)	\$13,290	\$9,967
SANYMPR-1412R-PE	1,359L with 10 roll-out drawers	\$14,908	\$11,180
SANY MPR-48B1-PW	Remote Alarm /Power Failure Alarm	\$382	FOC

- Double self closing doors (with viewing windows) and total capacity: >1,355 litres
- Automatic Cyclic defrost (as required) with no temperature spikes
- Temperature range: adjustable from 2°C to 23°C
- Alarms for high temperature, low temperature and door open
- Fast temperature recovery from door openings
- Large easy-to-read LED display of actual temperature (in 0.1°C increments)
- 5-Year Warranty



SANYMPR-1412R



SANYMPR-1412

VALID UNTIL:
31 May 2023
Ref. Code:
SANY-1412

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

