

Launch Promotion Microtube Shaker

Compact Benchtop Shaking Incubator

Code	Description	Price
AXYG I-4011	Axygen Microtube Shaker LCD -15°C to 100°C 230V 50/60 Hz	\$3,999

NEW



The Axygen Microtube Shaker allows uniform heating, cooling, and mixing of precious life science samples with a touch of the intuitive control panel.

Aluminum blocks-compatible with PCR plates/centrifuge tubes/microcentrifuge tubes are available for every life sciences applications.

The blocks utilize a magnetic adhesion technology and are self-recognized once exchanged.

- Precise temperature control for heating and cooling: Temp. range at -15°C ambient to 100°C
- Calibration mode ensuring accurate temperatures within $\pm 0.5^\circ\text{C}$ after every block exchange

FREE Extras

Purchase the Microtube Shaker and receive 1 of the below packages FREE:

Package 1 Microtube:

Choose either:



Heating Block,
cone bottom
holds 24 x
1.5ml tubes
AXYGI-4010-B

Heating Block,
round bottom
holds 24 x 2ml
tubes
AXYGI-4010-C

plus



AXYGMCT-150-C

2 x Cases 500 tubes per pack, 10 packs per case 1.5ml MaxyClear snaplock microcentrifuge tubes

Package 2 PCR Strip Tubes:

Choose 2 boxes:



Heating Block,
0.2ml tubes or
96-well PCR
microplates
AXYGI-4010-H

plus



AXYGPCR-0208-C
AXYGPCR-0208-CP-C
AXYGPCR-02-C
AXYGPCR-02CP-C
AXYGPCR-02-FCP-C
AXYGPCR-05-C
AXYGPCR-2CP-RT-C
AXYGPCR-96M2-HS-C



VALID UNTIL: 30 Sep Ref. Code: AXYG-NEW

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Launch Promotion Microtube Shaker




Heating Blocks

available at an additional cost



PCR Assay Tubes/Plates

Block	Description	Speed (rpm)	Well Dimensions (diameter x depth)
 AXYGI-4010-H	0.2ml tubes or 96-well PCR microplates	1,500	-
 AXYGI-4010-I	96 or 384-well microplates	1,500	-
 AXYGI-4010-D	8 x 5ml tubes	1,500	13.6 x 42mm
 AXYGI-4010-E	Cone bottom 8 x 5ml tubes	1,400	17 x 51.5mm

Microcentrifuge Tubes

Block	Description	Speed (rpm)	Well Dimensions (diameter x depth)
 AXYGI-4010-A	24 x 0.5ml tubes	1,500	8.1 x 26mm
 AXYGI-4010-B	Cone bottom 24 x 1.5ml tubes	1,500	10.9 x 30.8mm
 AXYGI-4010-C	Round bottom 24 x 2ml tubes	1,500	10.9 x 30.8mm

Centrifuge Tubes

 AXYGI-4010-F	8 x 15ml tubes	800	16.9 x 102.24mm
 AXYGI-4010-G	4 x 50ml tubes	600	29.6 x 105.95mm



VALID UNTIL: 30 Sep Ref. Code: AXGY-NEW

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.