

DURAN® YOUTILITY® Bottles

A New Generation Of Glass Laboratory Bottles

Code	Size	Price	Special
Clear YOUTILITY® Bottles with cyan cap and ring Pack of 4			
SCOT218812854	125ml	\$51.90	\$45
SCOT218813653	250ml	\$58.40	\$49
SCOT218814452	500ml	\$65.20	\$55
SCOT218815457	1Litre	\$90.80	\$79
Amber YOUTILITY® Bottles with cyan cap and ring Pack of 4			
SCOT218862859	125ml	\$251.10	\$215
SCOT218863658	250ml	\$300.30	\$249
SCOT218864457	500ml	\$320.50	\$275
SCOT218865453	1Litre	\$451.60	\$379



Pack of mixed colour bottle tags also available

The specially shaped gripping zones on both sides of the bottle permit easier and safer handling. With the new glass thread, the opening or closing of the YOUTILITY® bottle is significantly faster. Yet the thread remains fully compatible with DIN GL 45 closures and other accessories.

The slimmer YOUTILITY® bottle shape allows a more optimal use of limited space in autoclaves and laboratory refrigerators. The bottles can be conveniently stacked horizontally on their sides. A pre-defined labeling area is compatible with the dedicated self-adhesive YOUTILITY® bottle labels.



VALID UNTIL:
31 Aug 2023
Ref. Code:
SCOT-YOU



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.