

Parafilm® M Laboratory Film

200% Stretchability & Non Toxic

Code	Description		Price	Special
MARI7403810	38M x 10cm	1 Roll	\$69	\$49
MARI7407505	75M x 5cm	1 Roll	\$79	\$59

Use Parafilm® M to seal beakers, flasks, weighing bottles or graduated cylinders. This highly elastic, stretchable laboratory film prevents contamination, evaporation and spillage from fallen flasks. It stretches up to 200 % its original size and clings around irregular shapes and surfaces.

Parafilm® M complies with the provisions of the US FDA Food and Drug Administration when used below 55 °C and in accordance with GMP (Good Manufacturing Practice).



Physical Properties | Chemical Resistance

VALID UNTIL:
31 Jul 2023
Ref. Code:
MARI-Film



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.