




Optimum Wash
Pressure



EcoDry via Patented
AutoOpen Function²



Optimum Cleaning
Performance

Laboratory Washer/Disinfector With AD Rinse Water Pump

Was \$16,016 **Now \$14,250**




Code: MIEL62858306

- Two wash levels configured for washing wide neck and narrow neck glassware
- Spray arm and spray pressure monitoring ensuring a clean wash
- Compact: 600mm x 600mm x 850mm
- Underbench or standalone model
- EcoDry passive drying
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Dispensing pumps for liquid detergent and liquid neutraliser
- 2-year warranty



Baskets and inserts sold separately.

Please contact us for assistance in configuring your complete washing solution.



VALID UNTIL:
30 Jun 2023
Ref. Code:
MIEL-PG8583



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.