



Sentinel™ Gold Microprocessor Control System

Biological Safety Cabinet Class II 1.2m

Energy-Efficient | Quiet | Compact

Code	Description	Price	Special
ESCOAC2-4E8-SGT	Airstream E Class II Cabinet; 1.2m Type A2 with single piece work surface	\$17,397	\$10,840

Includes **FREE** adjustable height stand (710-860mm) with lockable castors, UV light and 2 power sockets (valued at \$2,154)



Certified by PHE (Public Health England)
for conformance to EN 12469

When safety is paramount!
The manual sash automatically activates internal lighting and airflow when set to the marked working height.

- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- Night Setback Mode reduces energy consumption by a further 60%
- Safety ensured : Interlocked UV light with timer
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Raised armrest external to cabinet for comfortable working position
- Recessed worksurface contains accidental spills
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 3-year warranty

VALID UNTIL:
30 Jun 2023
Ref. Code:
ESCO-AC2-4E8



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.