

Twister Magnetic Stirrer

The Tiny Hero!

Magnetic Stirrer 5L capacity, auto reverse/intermittent rotation, 100 - 3000rpm, vessel detection auto start

Code	IKAA20027010	IKAA20027388	IKAA20027389	IKAA20027390
No. of Driver Units	1	2	3	4
Price	Was \$1,070	Was \$2,010	Was \$2,945	Was \$4,020
Special	Now \$810	Now \$1,525	Now \$2,235	Now \$3,050

Is it a magnetic stirrer with up to 5L capacity, 3000 rpm with auto reverse and intermittent rotation?

Is it a single multi place stirrer for 9 tubes (Ø 23 or 28mm tubes) or a modular system for up to 30 stirrers?

Is it a shaker with adaptability for single vortex or for plates and multiple test tubes?

Yes. Yes. Yes. And more!

Robust with high protection class (IP66) and glass surface for corrosion resistance.

Auto detection of vessels for immediate stirring with flea.

This tiny hero will be the star of any lab!



[watch >](#)

Attachments

Code	Description	Price	Special
IKAA20103426	TW.IX Multi-position stirring attachment, 9 stirring positions for vessels with volume 8 - 60ml. Stirring position spacing: 30mm Max. ø vessels: 27.5mm	\$705	\$570
IKAA20110145	TW.VX Shaker / vortex attachment, Speed: up to 3000rpm, Shaking stroke: 4.5mm	\$565	\$430



IKAA20103426



IKAA20110145

VALID UNTIL:
30 Jun 2023
Ref. Code:
IKAA-Twist

Accessories also available.
Talk to us to find the right one to suit your needs



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.