

Laboratory Underbench Refrigerators and Freezers

Code	Description	Price	Special
ROLLLSFSR158	LSFSR158 Refrigerator, Solid Door, 158L	\$3,600	\$3,325
ROLLLSFSF98	LSFSF98 Freezer, Solid Door, 98L	\$3,980	\$3,660



LEC Underbench Laboratory Refrigerators and Freezers have been specially designed for storing combustible supplies safely.

- Compact footprint for space-constrained labs
- Dimensions 595 x 660 x 850mm (W x D x H)
- Freezer Temp control range: -18°C to -25°C
- Fridge Temp control range: +2°C to +10°C
- 12-month warranty as standard
(Extended 3-year warranty available)



Spark-Free



Access Port



Anti-Microbial Handle

VALID UNTIL:
30 Jun 2023
Ref. Code:
ROLL-LSFS



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.