

400L Laboratory Refrigerator

Stable Temperature



Code	Description	Price	Special
SANYLPR-400-PE	400L Upright with 5 adjustable height shelves	\$3,210	\$2,568



Designed for general laboratory storage with user-friendly features

The combination of temperature control accuracy, interior temperature uniformity, quick recovery and resistance to high ambient temperature delivers a quality storage refrigerator.

- Forced air circulation ensures improved temperature uniformity
- Two layer glass window in door for easy sample viewing
- Low noise design; 48dB(A)
- Access port on rear wall of cabinet for remote probe/cable access
- Outer door has cylinder keylock for security
- Cyclic Defrost (timer controlled)
- Easy to read display with digital microprocessor control
- Suitable for operation on 230V 50Hz 10 amp supply
- 3-Year Warranty



Visual & audible alarms notify end user of temperature deviation outside "safe" range

VALID UNTIL:
31 Jul 2023
Ref. Code:
SANY-LPR400



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.