



TwinGuard ULT 729L Freezer

Code	Description	Price	Special
SANYMDF-DU702VX-PE	Twinguard 729L Upright Freezer with 3 adjustable height shelves	\$26,500	\$19,975

-86°C Upright Freezer with two independent refrigeration circuits making this one of the most reliable ULT freezers on the market today!

If one refrigeration circuit fails – the other circuit maintains -70°C

- Small Footprint with VIP Plus vacuum insulation panel technology
- “Chamfered Design” ensures reduced installation space allowing more units to be placed side by side
- Storage for 576 x 2” cryo boxes with optional freezer racks
- Power failure, door, high and low temperature alarms for protection of valuable product
- Hot gas bypass around door frame and 5-fold gasket greatly reduce frost/ice build-up
- Automatic vacuum release port allows re-entry to freezer shortly after previous door closure
- Large colour touch panel display for clear indication of freezer performance
- USB port allows logged data of freezer’s operational status to be easily transferred to PC
- Outer door key lock and display lock for security
- 5-year warranty



Related Products

Cryo boxes and racks in a range of sizes available from various suppliers

VALID UNTIL:
30 Jun 2023
Ref. Code:
SANY-ULT729L



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.