

Single and Multi Channel Pipettors And Universal Tips

Axygen® Axypet® Single- and Multi-channel Pipettors

High levels of accuracy and precision with design innovations to make pipetting more comfortable during prolonged periods of time.

Description	Code	Price	Special
Single Channel			
0,5-10µL	AXYGAP-10	\$426	\$220
2-20µL	AXYGAP-20		
5-50µL	AXYGAP-50	\$384	\$230
100-1000µL	AXYGAP-1000		
Multi Channel			
20-200µL 8-Channel	AXYGAP-8-200	\$960	\$550



Axygen® Axypet® Pro Single- and Multi-channel Pipettors

Have been engineered to provide the highest levels of ergonomics, accuracy, and precision. All pipettors feature smooth plunger movement and extremely low pipetting forces to reduce the wrist strain and fatigue (RSI). Colored pushbuttons (included) can help define user or application for sample safety and lower risk of cross contamination.

Description	Code	Price	Special
Single Channel			
0.1-2µL	AXYGAP-2-P		
0.5-10µL	AXYGAP-10-P		
100-1000µL	AXYGAP-1000-P	\$357	\$250
5-50µL	AXYGAP-50-P		
10-100µL	AXYGAP-100-P		
Multi Channel			
0.5-10µL 8-Channel	AXYGAP-8-10-P	\$969	\$515



Universal Tips are designed to fit a wide variety of single and multi-channel pipettors, compatible with most popular brands of pipettors.

These polypropylene tips are nonpyrogenic and RNase-/DNase-free.

Our Axygen tips are available in Extended Length, Maxymum Recovery (low retention), Eppendorf-style, Wide-bore, Sterile and Non-sterile, Filtered and Non-filtered, Racked, Refillable Rack System and Bulk to suit your needs.

Flexible walls and a series of internal sealing rings ensure a secure fit with less force required to load and eject.

Description		Code	Price	Special
Filter Racked				
50µL Sterile	4800 items	AXYGTF-50-R-S	\$710	\$419.95
100µL Universal Fit Sterile	10x96 items	AXYGTF-100-R-S	\$138	\$66.95
200µL Max Rec P200 Rack Sterile	10X96 items	AXYGTF-200-L-R-S	\$124	\$88.95
Universal Racked				
200µL Yellow	10x96 items	AXYGT-200-Y-R-S	\$74	\$43.95
Universal Bulk				
200µL Yellow	1,000 items	AXYGT-200-Y	\$25	\$16.95
1000µL Blue	1,000 items	AXYGT-1000-B	\$30	\$18.95
Maxymum Recovery Range				
10µL Reach MR Rack Sterile	960 items	AXYGTXL-10-L-R-S	\$131	\$84.95
1000µL Clear Universal Fit Bulk	1,000 items	AXYGT-1000-C-L	\$42	\$28.95



VALID UNTIL:
31 Aug 2023
Ref. Code:
AXYG-Pet/Tips



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.