



Spray Arm  
Sensing



Optimum Cleaning  
Performance



Crevice-free  
Chamber

## PG8583 CD Washer/Disinfector

Laboratory Washer/Disinfector with DryPlus HEPA filtered drying.


- 14 standard programs + 2 freely configurable programs
- Two wash levels configured for washing wide neck and narrow neck glassware
- Conductivity monitoring to ensure load is clean
- Pull out drawer for housing detergent and neutraliser
- Spray arm and spray pressure monitoring ensuring a clean wash
- Compact: 900mm x 700mm x 835mm W x D x H (820mm with lid removed)
- Underbench or standalone model
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Stainless steel interior/exterior for durability and cleanliness
- ADP model including pump for AD water and waterproofing kit
- Many options and accessories to select from
- Dispensing pumps for liquid detergent and liquid neutraliser
- 2-year warranty



Was \$29,439 **Now \$26,495**

Code: MIEL62858351

Baskets and inserts sold separately.  
Please contact us for assistance in configuring your complete washing solution.



VALID UNTIL: 31 Aug 2023

Ref. Code: MIEL-PG8583



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.