

Moisture Analysis Made Simple!

Ohaus Moisture Analysers consistently deliver accurate and reliable results

- Flexible Performance, for a wide variety of applications
- Simple to Operate, for faster, smarter and more efficient analysis
- Compact, Easy-to-Clean Design

Code	Description Moisture Analyser	Price	Special
OHAU30246943	MB120 Capacity: 120g; Readability: 0.01%/0.001g	\$5,933	\$4,450
OHAU30246944	MB90 Capacity: 90g; Readability: 0.01%/0.001g	\$4,933	\$3,700
OHAU80252472	MB25 Capacity: 110g; Readability: 0.05%/0.005g	\$3,322	\$2,490

MB120

Delivers functionality and precision that you can rely on. Offering advanced features that simplify processes whilst significantly improving lab throughput/efficiency, and robust housing, which is attractive, durable and easy to clean.

Key Features:

- SmartGuide™ software: a revolutionary feature that automatically analyses a sample and creates a method, thereby removing the guesswork from the initial setup
- Temperature Range: 40 to 230°C using halogen lamp
- 4 drying profiles and 7 shut-off criteria
- 110mm TFT touchscreen
- Internal memory capable of storing 100 methods and 1000 measurements



VALID UNTIL: 31 Aug 2023
Ref. Code: OHAU-MB

Moisture Analysers

MB90

Featuring a precisely-controlled halogen heating system which helps to yield fast and repeatable results with readability up to 0.01%, with an intuitive touchscreen icon-driven menu guides users through every step of the operation process.

Key Features:

- Temperature Range: 40 to 200°C using halogen lamp
- 2 drying profiles and 3 shut-off criteria
- 110mm TFT touchscreen
- Internal memory capable of storing 2 methods



MB25

An ideal combination of results and value, the MN25 is extremely user-friendly. Its compact and easy-to-clean design allows it to fit seamlessly into any workspace. Great for routine measurements!

Key Features:

- Temperature Range: 50 to 160°C using halogen lamp
- 1 drying profile and 3 shut-off criteria
- LCD Display



Related Products

Accessories / Consumables:

Code	Description
OHAU80252478	Pans, 3 Reusable, 7mm
OHAU80252479	Pans, 3 Reusable, 14mm
OHAU30585411	Pans, 50, Disposable, Aluminium
OHAU80850087	Filter, 200, Glass Fiber Pads

VALID UNTIL:
31 Aug 2023
Ref. Code:
OHAU-MB



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.