

Aquastab Algaecide from Scharlau

**Reduce water bath contamination with
Aquastab Algaecide!**



Marketed by Scharlau

Was \$115.00

Now \$89.95

Code

SCAR0000LZB929

Aquastab Algaecide is a highly effective blue, clear additive to reduce the contamination from bacteria and fungi in your water bath.

Use 5ml/10L Water - 100ml bottle makes up 200L

Highly effective algicide only to be used with water in bath thermostats.

Composition:

2% Poly(dimethyl imino)-2-hydroxy-1,3-propanediyl chloride

Temp. Range: 5°C to 90°C

[Download](#) Safety Data Sheet

VALID UNTIL:
31 Aug 2023

Ref. Code:
SCAR_Aqua



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.