

Pioneer PX224 Analytical Balance

**Affordably Priced
Intuitively Designed**



Was \$3,000

Now \$2,250

Code: OHAU30428792

The Ohaus Pioneer PX224 offers competitive performance, with high accuracy and reliable results – at a great price!

Features:

- Maximum Capacity: 220g
- Readability: 0.1mg
- Pan Size: 90mm
- Built-in static removal bar
- 3-year warranty

Design Features

InCal™ - Semi automatic internal calibration, user-selectable environmental filters and brightness settings, auto-tare, RS232 and USB ports included with user selectable communication settings and data print options, security bracket, GLP/GMP data output with real-time clock, stability indicator, four operating languages.

Need communication?
Add an SF40A dot-matrix
printer for an additional
\$640.00!



Code: OHAU30045641

VALID UNTIL:
31 Aug 2023
Ref. Code:
OHAU-PX224



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.