

165L IncuSafe CO₂ Incubator

Was \$10,796


Now \$8,595

Code: MCO-170AICUVL-PE

3 In Stock

CO₂ Incubator with IR sensor and UV decontamination

- Corrosion-resistant stainless steel copper enriched alloy interior (inCu-saFe)
- Inner glass door for viewing cultures
- Temperature: Ambient +5°C to 50°C (±0.1°C)
- CO₂ : 0 to 20% (±0.15%)
- Direct heat and air jacket
- 4 adjustable shelves
- PID control with colour LCD display
- 5-year warranty



VALID UNTIL:

30 Sep 2023

Ref. Code:

SANY-CO2170L

Price available for current stock only



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.