

**VIP ECO -86°C  
ULT Upright  
528L Freezer**

**Preservation**  
sub-zero

Was \$22,839

**Now \$17,975**

Code: MDF-DU502VH-PE

This upright freezer is utilizing ECO "Green Gas" refrigerants combined with inverter compressors for leading class performance, low energy consumption and low running costs.

- Storage for 384 x 2" cryoboxes with optional freezer racks
- Power failure, door, high and low temperature alarms for protection of valuable product
- Low noise design; <52dB(A)
- Dry remote alarm contacts and access ports
- Suitable for operation on 240V 50Hz 10 amp supply
- Insulated inner doors with eagle latches keep the cold in
- 1,993mm tall – fits under 2m door frame
- Hot gas bypass around door frame and 5-fold gasket greatly reduce frost/ice build-up
- Automatic vacuum release port allows re-entry to freezer shortly after previous door closure
- "Chamfered Design" ensures reduced installation space allowing more units to be placed side by side
- Large colour touch panel display for clear indication of freezer performance
- USB port allows logged data of freezer's operational status to be easily transferred to PC
- Outer door key lock and display lock for security
- 5-year warranty



**watch >**

**Related Products**

Cryo boxes and racks in a range of sizes available from various suppliers

VALID UNTIL:

30 Sep 2023

Ref. Code:

SANY-DU502VH



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in NZD excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.