


## Isotherm Laboratory Benchtop Incubators

Code	Size	Price	Special
IFA-110-8	110L	\$5,210	<b>\$3,995</b>
IFA-170-8	170L	\$6,847	<b>\$4,995</b>
IFA-240-8	240L	\$8,463	<b>\$6,500</b>

These forced convection incubators are used for high-volume thermal convection applications.

They provide uniform temperature throughout the chamber necessary for annealing, drying, sterilizing, and other lab functions.



### Pre-heat Chamber Technology

Ensures stable heating and maximum temperature uniformity in the chamber



### Easy to Clean

- Ambient + 7.5°C to 100°C
- Dual display showing set and actual temperature
- Inner glass door for sample observation without affecting chamber temperature
- Isocide™ Antimicrobial Powder Coating: Eliminates 99.9% of surface bacteria within 24 hours of exposure
- Single-piece stainless steel chamber with directly mounted shelves making oven easy to clean
- 4-zone heated air jacket for maximum temperature uniformity
- Supplied with 2 shelves (additional shelves can be ordered)
- Permanently lubricated maintenance-free fan
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10 program memory



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
ESCO-IFA

VALID UNTIL:  
30 Nov 2023