

Low-Speed Centrifuge Clinical Package

Reliable & Compact

Only \$4,995

Usually \$5,842

Code: GYROGZ-0416-CLINICAL-PK



A range of alternate rotors/adaptors are available,
contact us for more information

Package Includes

- 1 x 416 compact low-speed centrifuge
- 1 x Swing out rotor, 4 x 100ml
- 4 x Round 100ml buckets with Biosafety Cap
- 4 x Adapters, 5 x 2.6-7ml tubes

This package will quickly prove to be a valuable asset to any clinical workflow.

Dependable performance, essential safety features, versatility, and ease of use – everything you need and more!

Features:

- Imbalance cut-off and concurrent alarms (imbalance, speed and over-heating) for safe operation
- Unique air flow design and anti-vibration structure to minimize heat and noise generation
- Soft Start/Stop function to minimize agitation (slow run-up/down)
- Max. Speed: 4,000rpm ; Max Force: 2,774 x g
- Highly configurable for application diversity
- Very quiet operation at less than 56 dB, allowing a comfortable laboratory environment
- Lid drop protection allows convenient and safe loading or unloading of samples
- Automatic lid opening when rotor is static, for easy retrieval of samples
- Up to 10 programs memory for storing routine protocols
- 2-year warranty



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
GRYO-416

VALID UNTIL:
15 Dec 2023