



Crevice-free  
Chamber



Optimum  
Cleaning



Optimum  
Wash Pressure



Short  
Cycle Times

## Professional Laboratory Washer

Are you still manually washing glassware?  
Make better use of your staff's valuable time.

Code	Description	Price	Sale
MIEL62850422	PG8504 ADP SS 230V STAINLESS STEEL, 60CM	\$11,749	<b>\$10,570</b>

Laboratory Underbench/Free-Standing Washer with stainless steel interior combining powerful 500L/min pump with integrated heating elements enabling safe washing cycles with temperatures up to 70°C.



Miele offers enhanced laboratory safety by removing the manual washing step from your wash process.

There is reduced likelihood of glass breakage that can result in cuts to lab staff.

### Miele - good for occupational health and safety

- Door dispensers for powder detergent and neutraliser
- Suitable for operation on 240V 50Hz 13 amp supply
- Underbench or standalone model
- 6 standard programs + 2 vacant program positions to meet your wash requirements
- Touch on steel controller with 3-line display for easy viewing of washer status
- Stainless steel interior/exterior for durability and cleanliness

Many options and accessories to select from,  
contact Bio-Strategy for assistance in configuring your washer.



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
MIEL-ADP

VALID UNTIL:  
30 Nov 2023