

82L Pharmaceutical Fridge

Practical & Fit-for-Purpose

Now \$2,560

Was \$2,890

Code: ROLLPEGR82
Glass Door

[View content at a glance](#)



Now \$2,300

Was \$2,630

Code: ROLLPESR82
Solid Door

[Protect light-sensitive samples](#)



The Lec Under/On Bench Pharmaceutical fridges are compact, reliable, and efficient. Perfect for personal or low volume storage of medical or laboratory samples.

Features include:

- 2 x locks as standard – offering security and controlled access
- Temperature Control: +2°C to +8°C
- Gross Capacity: 82 litres
- External Dimensions (WxDxH): 495 x 515 x 655mm



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
ROLL-82LFridge

VALID UNTIL:
15 Dec 2023