


# CMT6 Coverslipper

**\$25,000**

**Normally \$48,300**

Code: EPRN970010

**Ex-Demo Unit**



## Efficiency for Time Savings

- Coverslip hundreds of slides per hour with a self-cleaning dispensing needle
- Quick and easy change of different sizes and styles of coverslips for versatility with different slides and tissue types

## User Safety

- Load and unload rack without opening the hood
- Internal vapor extraction system with fan and permanent air-flow avoids any leakage of vapors at openings
- Integrated charcoal filtration
- Tool-free cleaning and maintenance

**Compact** - Designed to fit almost any lab space and under a 77cm or larger hood

## Compatible with:

### Coverslips

EPBRCS24401GP	24mm x 40mm
EPBRCS24501GP	24mm x 50mm
EPBRCS24601GP	24mm x 60mm

**Discounts available - Contact us for pricing**



Ref. Code:  
EPRN-CTM6

VALID UNTIL:  
31 Mar 2024  
or until sold

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.