

# 1.2m Airstream Horizontal Laminar Flow Cabinet



## Energy efficient DC ECM Motor

(70% energy savings  
compared to AC motor)

One of the quietest (only 55.8dB(A)) and safest cabinet choices  
for sample protection requirements.

**Only \$5,795**

Was \$7,671

### Includes:

*for current stock only*

Code	Description
ESCOLHG-4AG-F8	Horizontal Laminar Flow Cabinet 1.2m
ESCOEO-HE	Electrical Outlet x 2
ESCOUV-30A	Safety Interlocked UV light with timer
ESCOFC-LHG/S-4A	Front Cover

**In Stock**



Optional adjustable height stand available for  
only **\$895** (usually \$1,515)  
just add code ESCOSTC-4AO

- Fixed height polycarbonate sash and UV resistant glass side panels allow entry of "natural" light.
- ULPA filter (10x filtration efficiency than of HEPA filter) delivers ISO Class 3 air cleanliness.
- Sentinal™ Gold microprocessor controller with digital display showing cabinet status.
- Stainless steel spill retaining worksurface for ease of cleaning
- Stable and self-compensating airflow.
- Real time airflow monitoring with alarms.
- Standby mode provides more energy savings and maintains an ISO Class 5 work zone.
- ISOCIDE™ Powder Coating – painted surfaces have silver ion impregnated coating to inhibit microbial growth.
- 3-year warranty



Ref. Code:  
ESCO-Flow

VALID UNTIL:  
31 Mar 2024

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.