

1.2m Airstream E Class II Microbiological Safety Cabinet

When safety is paramount!
Operator protection, product protection and environment protection!

Only \$11,325

Was \$15,673

Includes:

for current stock only

Code	Description
ESCOAC2-4E8	Class II Microbiological Safety Cabinet, 1.2m
ESCOEO-HE	Electrical Outlet x 2
ESCOUV-30A	Safety Interlocked UV light with timer

In Stock



Optional adjustable height stand available for
only **\$895** (usually \$1,515)
just add code ESCOSTC-4AO

- UV resistant glass side panels allow entry of "natural" light
- Manual sash automatically activates internal lighting and airflow when set to the marked working height
- Night set back mode reduces power consumption by a further 60%
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness.
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Raised armrest external to cabinet for comfortable working position
- Multi-piece stainless steel work surface is easy to remove for cleaning
- Certified by PHE (Public Health England) for conformance to EN 12469
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coating – painted surfaces have silver ion impregnated coating to inhibit microbial growth.
- 3-year warranty



Ref. Code:
ESCO-E-ClassII

VALID UNTIL:
31 Mar 2024

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.