

## PG8583 Washer/Disinfector

**Only \$14,365**

Was \$15,965

Code: MIEL62858306  
*accessories not included*

- Underbench or stand-alone model
- Stainless steel interior/exterior for durability and cleanliness
- Spray arm and spray pressure monitoring ensuring a clean wash
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Passive EcoDry drying function – glassware can be dried
- ADP model including pump for AD water and waterproofing kit
- Two wash levels
- Many options and accessories to select from
- Dispensing pumps for liquid detergent and liquid neutraliser
- 2-year warranty



**15% discount** on optional accessories  
Contact us to discuss your complete washing solution

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
MIEL-PG8583

VALID UNTIL:  
31 Mar 2024