

Model B150-BASE

150L Top Loading Autoclave

In Stock
1 Unit Only!

Only \$28,800

Was \$36,035

Code: PRIOEV150-BASE

Includes:

3 baskets &
Export Crate



Complying with AS/NZS standards, the unit has grade 316L stainless steel chamber and is inherently safe with thermal and pressure interlocks.

- Single handed lever action door mechanism
- 5 program memory
- Water cooled drain protection as standard
- Quickseal single-action door fitted with thermal and pressure locks to prevent opening at load temperature above 80 °C and pressures above 0.2 bar
- TACTROL 3 microprocessor control system with simple and fully variable setting of process time, temperature and with graphic indication of cycle status
- Operation up to 140C - 2.5 Bar
- Large 150L cylindrical chamber within a compact frame
- Forced air cooling as standard helping to reduce cycle times
- Automatic timed free steaming as standard for improved air removal and sterilising performance
- Steam generated within the chamber (electrically heated model)
- Complete with thermocouple entry and test ports
- Delayed start and media warming function as standard for ready to pour media at the start of the day
- Low water level sensor
- Epoxy coated panels and frame members are treated with an anti-bacterial agent which is effective against all bacteria and fungi including MRSA
- Pressure vessel carries a 10-year warranty and is insurance approved

Hurry won't last long!

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
PRIO-EV150

VALID UNTIL:
31 Mar 2024
unless sold prior