



## -40°C Biomedical Eco Freezer with Natural Refrigerant 479L



**Only \$7,600**

Was \$9,600

Code: MDF-MU549DH-PE

With a **free** set of storage bins  
(value \$1,446.00)

*price valid for current stock only*

**2 Units  
in Stock**



- Capacity: 479L (2 x 240L)
- Quiet in operation : 42dB(A)
- Includes height adjustable polyethylene coated steel shelves
- Push button manual defrost
- Alarm test button
- Two access ports as standard (1 per chamber)
- Castors for mobility and locating feet when in working position
- Keyed door locks for security
- High temperature, low temperature and power failure alarms
- Dry remote alarm contacts
- Inverter compressor technology with variable speed compressors ensure the most efficient energy use and reduced heat output
- Suitable for operation on 220/240V 50Hz 10 amp supply
- Drain hose to remove defrosted water
- 5-year warranty

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
SANY-ECO-479L

VALID UNTIL:  
31 Mar 2024