

# Isotherm Laboratory Benchtop Incubators

Fan Forced, Ambient + 7.5°C to 100°C

Code	Size	Price	Special
IFA-54-8	54L	\$3,976	<b>\$3,100</b>
IFA-170-8	170L	\$6,847	<b>\$4,480</b>
IFA-240-8	240L	\$8,463	<b>\$5,540</b>


Supplied with 2 shelves  
(additional shelves can be ordered)

**In Stock**



## Designed For Laboratory Culturing Applications

- Isocide™ Antimicrobial Powder Coating: Eliminates 99.9% of surface bacteria within 24 hours of exposure
- Single-piece stainless steel chamber with directly mounted shelves making incubator easy to clean
- Permanently lubricated maintenance-free fan
- Air flow slider to adjust exchange rate of air
- Microprocessor PID control with 10 program memory



## Pre-Heat Chamber Technology

4-zone heated air jacket for maximum temperature uniformity

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: ESCO-ISO\_IFA  
VALID UNTIL: 30 Apr 2024  
*price valid for in stock items only*