



- Equipped with energy-efficient ECM blower from ebm-papst Germany
- Extremely low energy consumption (190 watts)
- Quietest cabinet in the industry at 51 dBA

NordicSafe Class II Microbiological Safety Cabinet

The Industry's
Most Comfortable
Energy Efficient Cabinet

Certified by PHE (Public Health England)
for conformance to EN 12469



Optional adjustable height stand for
\$995 (usually \$1,425)
just add code ESCOSTC-4CO

Only \$14,000

Was \$19,590

Code: ESCONC2-4L8

Included **FREE**
UV light & 2 power sockets

- Motorised sash ensures sash travels to working height automatically
- Ergonomic sloping front combined with UV resistant tempered glass side panels allow entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Audible and visual alarms
- Raised armrest external to cabinet for comfortable working position
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- Night set back mode reduces power consumption by a further 60%
- 4-year warranty

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: **ESCO-NC2-4L8** VALID UNTIL:
30 Apr 2024

price valid for in stock items only