

# SensiStrip LFD Kits (Lateral-Flow-Device)

**Introductory Offer**  
**Just \$430 per kit**

**Save time,  
without compromise**

The *SENSiStrip* Lateral Flow Device (LFD) Allergen kits are the perfect tool to control any potential cross- contamination with food allergens, avoiding external analyses and thus considerably reducing the time needed to release a food product batch to the market.



- Quick & Easy to use
- Highly-specific immunoassay-based kits
- Sensitive, with low levels of detection
- Validated for food matrices, rinse water and equipment



Sample preparation


Transfer to incubation vial

Place strip into vial


Evaluate strip



3 min



5 min



Ideal for use with the LFD Reader  
See current **[promotion](#)**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
GOLD-Sensi

VALID UNTIL:  
31 Mar 2024

These rapid tests, with high sensitivity and specificity, include all the necessary material to perform 20 individual analyses per kit, and are validated on food matrices, swabbing and rinse water.

Also available in a 5-strip pack size.

### SENSIS*Strip* LDF Kits 20 Pack

Code	Description	Code	Description
GOLDHU0030110	Almond	GOLDHU0030129	Pecan Nut
GOLDHU0030112	Casein	GOLDHU0030130	Pistachio
GOLDHU0030113	Cashew	GOLDHU0030131	Sesame
GOLDHU0030114	Coconut	GOLDHU0030132	Total Soy
GOLDHU0030115	Crustaceans (Tropomyosin)	GOLDHU0030133	Walnut
		GOLDHU0030111	Beta-Lactoglobulin
GOLDHU0030116	Egg White	GOLDHU0030117	Fish LFD
GOLDHU0030119	Hazelnut	GOLDHU0030120	Lupin
GOLDHU0030122	Macadamia Nut	GOLDHU0030121	Lysozyme
GOLDHU0030123	Milk Protein	GOLDHU0030126	Ovalbumin
GOLDHU0030125	Mustard	GOLDHU0030127	Brazil Nut
GOLDHU0030128	Peanut	GOLDHU0030118	Gluten

Swabbing kit also available for customers who conduct surface testing – please ask for details & pricing.



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
GOLD-Sensi

VALID UNTIL:  
31 Mar 2024