

# Hub RN Needles

Recently purchased a Hamilton syringe with removable needle and looking for spares or want to re-stock?

**Check out our special pricing across all standard small and large hub RN needles**

**From \$81**

Removable needles give you the best of both worlds

1. the benefit of a replaceable needle for the occasional bend or clog
2. the minimal dead volume of a cemented needle. The unique removable needle design installs the needle precisely at the zero line

**Large**



Compatible with syringes ranging from 250uL up to 10mL.

**Small**



Compatible with syringes ranging from 2.5uL up to 100uL.

**Need custom needles, just give us a call**

*see pricing list page 2*

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
HAMC-HubRN

VALID UNTIL:  
30 Sep 2024

## Hub RN Needles

### Packs of 6

Code	Size	Price	Special
<b>Large</b>			
HAMC7779-02	26s/2"/2	116	<b>87</b>
HAMC7779-04	26/2"/2		
HAMC7779-03	22s/2"/2	156	<b>117</b>
HAMC7779-01	22/2"/2	116	<b>87</b>
HAMC7780-01	26s/2"/3	108	<b>81</b>
HAMC7780-02	26/2"/3		
HAMC7780-03	22s/2"/3	163	<b>122</b>
HAMC7780-04	22/2"/3	119	<b>89</b>
HAMC7784-03	26s/2"/5	248	<b>186</b>
HAMC7784-04	26/2"/5		
HAMC7784-01	22s/2"/5	288	<b>216</b>
HAMC7784-02	22/2"/5	248	<b>186</b>
<b>Small</b>			
HAMC7758-02	26s/2"/2	128	<b>96</b>
HAMC7758-04	26/2"/2		
HAMC7758-03	22s/2"/2	172	<b>129</b>
HAMC7758-01	22/2"/2	128	<b>96</b>
HAMC7762-06	33/1.5"/3	184	<b>138</b>
HAMC7762-05	32/2"/3		
HAMC7762-04	31/2"/3		
HAMC7762-03	30/2"/3	148	<b>111</b>
HAMC7762-02	28/2"/3		
HAMC7762-01	7/2"/3		
HAMC7768-01	26s/2"/3	108	<b>81</b>
HAMC7768-02	26/2"/3		
HAMC7770-01	22s/2"/3	148	<b>111</b>
HAMC7770-02	22/2"/3	108	<b>81</b>
HAMC7784-07	26s/2"/5	248	<b>186</b>
HAMC7784-08	26/2"/5		
HAMC7784-05	22s/2"/5	288	<b>216</b>
HAMC7784-06	22/2"/5	248	<b>186</b>



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
HAMC-HubRN

VALID UNTIL:  
30 Sep 2024