

Simple, Accurate & Reliable Moisture Analysis

Only \$4,450

Was \$5,933

Code: OHAU30246943

The MB120 makes the moisture determination process faster, easier, and more efficient.



watch >

Offering advanced features that simplify processes whilst significantly improving lab throughput/efficiency, and robust housing, which is attractive, durable and easy to clean.

- Simple to Operate, with a compact, easy-to-clean design
- Smarter features & SmartGuide™ software deliver faster and more efficient analysis
- Flexible Performance, for a wide variety of applications
- SmartGuide™ software: a revolutionary feature that automatically analyses a sample and creates a method, thereby removing the guesswork from the initial setup
- Capacity: 120g; Readability: 0.01%/0.001g
- Temperature Range: 40 to 230°C using halogen lamp
- 4 drying profiles and 7 shut-off criteria
- 110mm TFT touchscreen
- Halogen heating technology
- Internal memory capable of storing 100 methods and 1000 measurements

Accessories / Consumables:

Contact us for pricing

Code	Description
OHAU80252478	3 Reusable pans 7mm
OHAU80252479	3 Reusable Pans 14mm
OHAU30585411	50 Disposable Aluminium Pans
OHAU80850087	200 Glass Fibre Filter Pads

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
OHAU-MB120

VALID UNTIL:
31 Mar 2024