



## 165L CO<sub>2</sub> Incubators

**Only \$6,100**

**Was \$8,700**

Code: SANYMCO-170ACL-PE



TCD Sensor  
3 Adjustable Shelves

**Only \$6,650**

**Was \$9,656**

Code: SANYMCO-170AICL-PE



Infrared (IR) Sensor  
4 Adjustable Shelves

**In Stock**

The latest IncuSafe CO<sub>2</sub> Incubators offer outstanding quality in performance, maximise cell culture productivity and provide optimum results and reproducibility.

Choice of TCD sensor or dual wavelength single beam IR CO<sub>2</sub> sensor with continuous auto zero calibration.

- Temperature: Ambient +5°C to 50°C (±0.1°C)
- CO<sub>2</sub>: 0 to 20% (±0.15%)
- inCu-saFe interior: provides the germicidal properties of copper and the corrosion resistance of stainless steel without tarnishing!
- Inner glass door for viewing cultures without changing the environment
- Direct heat and air jacket for environment stability
- Humidity: 95%RH (±5%RH)
- Field reversible door
- Peltier powered dew stick to draw condensation away from inner glass door
- PID control system with colour touchscreen display
- USB port for data transfer
- 5-year warranty

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
SANY-MCO-170

VALID UNTIL:  
31 Mar 2024