

## Airstream E Class II Cabinet Package

**Only \$14,995**

Was \$21,892

Code: ESCOAC2-5E8

Includes: UV light and 2 power sockets

*Pricing valid for current stock only*

**1 Unit Only  
At This Super  
Low Price**

When safety is paramount!

Operator protection, product protection and  
environment protection!



Optional adjustable height stand available for  
only **\$995** (usually \$1,763)  
just add code ESCOSTC-5AO

- Quiet : only 53.6dB(A)
- manual sash automatically activates internal lighting and airflow when set to the marked working height.
- UV resistant glass side panels allow entry of "natural" light
- Ergonomic sloping front combined with tempered glass side panels allows entry of more natural light.
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%) improves safety and delivers ISO Class 3 air cleanliness.
- Sentinel Gold microprocessor controller with digital display shows cabinet status.
- Raised armrest external to cabinet for comfortable working position.
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient.
- Removable easy-to-clean paper catch.
- Energy efficient DC ECM Motor ( with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth.
- 3-year warranty

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:  
ESCO-AC2-5E8

VALID UNTIL:  
30 Apr 2024