

GSD Absorbance 96

The most compact
ELISA reader on the market!

Only \$8,500

Was \$10,780

Code: GOLDEIAREAD002



Robust performance and accurate
results - with a small footprint.

A new category of plate reader with the mission
to simplify the workflow in the laboratory.

Key Features:

- Compact
- Cost-effective
- Easy to use (plug-and-play, intuitive multi-language interface)
- 96 individual detection units enabling simultaneous signal detection
- Compatible to GSD plates and other brands
- No license fee for software use

Specifications

- Detection method: Absorbance
- Detection mode: Endpoint, kinetic
- Microplate types: 96-well microplates
- Wavelength selection:
Example combinations (FWHM 10 nm): 450, 492, 562, 620 nm / 405, 450, 540, 630 nm / 492, 562, 605, 650 nm. Other filter combinations in the range of 400–1000 nm
- Power supply: USB connection 5 V

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
GOLD-ELISA

VALID UNTIL:
31 May 2024