

# Rotary Microtome HM 325

## Section Safely With Excellent Quality

### Designed with the end user in mind

- High stability and precision for premium-quality paraffin sections
- Unique section thickness setting with removable knobs allowing operation with both the right and left hands
- Resettable electronic section counter
- Two quick-trim stages (10 µm and 30 µm) independent from the preset fine section thickness
- X/Y fine orientation with reproducible zero positioning
- Manual coarse feed wheel
- Specimen retraction
- Long stroke length (64 mm) allows high quality sections of Macro / SuperMega cassettes
- Large, removable section waste tray covers the entire working area

Code:	Description - Microtome HM325 including Universal Cassette Clamp
EPRN902100ER	with ER Blade Holder
EPRN902100A	with E Blade Holder



### Compatible with:

EPBRSF41296SP	Superfrost Plus Slides
EPBRSF41296	Superfrost Slides
EPBRSF41296PL	Superfrost Polysine
EPKM3053835	Ultra Microtome Blades

**In Stock**

**Request a quote today!**

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: EPRN-HM325  
VALID UNTIL: 30 June 2024