



- Equipped with energy-efficient ECM blower from ebm-papst Germany
- Extremely low energy consumption (190 watts)
- One of the quietest cabinet in the industry at 51 dBA

NordicSafe Class II Microbiological Safety Cabinet 1.2m

**Available for
Immediate Delivery**

Only \$14,415

Usually \$19,590

Code: ESCONC2-4L8

Includes: UV light, 2 power sockets
and adjustable height stand



Certified by PHE (Public Health England)
for conformance to EN 12469

- Motorised sash automatically travels to correct height and activates internal lighting and airflow
- Night set back mode reduces power consumption by a further 60%
- Ergonomic sloping front combined with tempered UV glass side panels allows entry of more natural light
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Audible and visual alarms
- Raised armrest external to cabinet for comfortable working position
- Multi-piece stainless steel worksurface is easy to remove for cleaning
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Removable easy-to-clean paper catch
- Energy efficient DC ECM Motor (with 70% energy savings compared to AC motor)
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 4-year warranty

pricing for current stock only

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: **VALID UNTIL:**
ESCO-NC2-4L8 **30 June 2024**