

1730R Refrigerated Microcentrifuge

Ex-Demo
1 Unit Only!

Just \$5,400

Usually \$8,292

Code: GYROGZ-1730R-24-PK

Includes

- 1 x Refrigerated Microcentrifuge 1730R
- 1 x biocontainment fixed angle rotor (24 x 1.5/2.0ml)

Excellent Condition
Comes with a 2-year warranty



Functional & Flexible

A "must-have" item for
molecular and cellular labs

Features:

- Max. speed 17,000 rpm (27,237 x g)
- Temp Control Range: -20°C to 40°C
- Fast cooling function to 4 °C in 5 minutes for temperature sensitive samples
- Very quiet operation, ≤56 dB
- "PULSE" key for quick spin
- Selectable time mode of "From Starting" or "At Set Speed"

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
GRYO-1730R_Demo

VALID UNTIL:
31 May 2024