

1700 Series Syringes

From \$114

Gastight syringes from Hamilton can be used to accurately dispense both liquids and gases.

Available in a range of termination styles, with working volumes from 1uL up to 500uL.

These specific 1700 series syringes have a removable hub termination and are compatible with large and small hub RN needles.

They have a precision-machined PTFE plunger tip which creates a leak-free seal.

With the tight fit, the tip essentially wipes the interior of the syringe barrel free of sample.

This feature is particularly useful with heterogeneous samples as it reduces the chance that a deposit will occur and cause the plunger to freeze.



see pricing list page 2

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
HAMC-1700_Series

VALID UNTIL:
31 May 2024

1700 Series Syringes

Syringe With Needle

Code	Model	Volumne	Size	Price	Special
HAMC80030	1701RN		(26s/2"/2)	\$150	\$120
HAMC80065	1701RNR	10ul	(22s/2"/3)		
HAMC80014	1701RN		(32/2"/3)	\$154	\$123
HAMC80230	1702RN		(22s/2"/2)	\$164	\$131
HAMC80265	702RNR	25ul	(22s/2"/3)	\$162	\$130
HAMC80930	1705RN		(22s/2"/2)	\$164	\$131
HAMC80965	1705RNR	50ul	(22s/2"/3)	\$162	\$130
HAMC81030	1710RN		(22s/2"/2)	\$166	\$133
HAMC81065	1710RNR	100ul	(22s/2"/3)	\$164	\$131
HAMC81130	1725RN		(22s/2"/2)	\$172	\$138
HAMC81165	1725RNR	250ul	(22/2"/3)	\$164	\$131
HAMC81230	1750RN		(22/2"/2)		
HAMC81265	1750RNR	500ul	(22/2"/3)	\$176	\$141

RN = Removable Needle | RNR = Removable Needle Rheodyne

Syringe Only - Needles Sold Separately

Code	Model	Volumne	Price	Special
HAMC7653-01	1701RN	10ul	\$143	\$114
HAMC7654-01	1702RN	25ul		
HAMC7655-01	1705RN	50ul	\$154	\$123
HAMC7656-01	1710RN	100ul		
HAMC7657-01	1725RN	250ul	\$160	\$128
HAMC7658-01	1750RN	500ul	\$179	\$139



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
HAMC-1700_Series

VALID UNTIL:
31 May 2024

0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com