

PG8593 Washer/Disinfector

Only \$21,927

Usually \$24,364


Code: MIEL62858306



- Underbench or stand-alone model
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- Active DryPlus HEPA filtered drying ensures your loads come out dry
- Two wash levels
- Many options and accessories to select from
- Dispensing pump for liquid detergent , optional dispensing pump required for liquid neutraliser
- 2-year warranty



Stainless steel interior/exterior
for durability and cleanliness



Spray arm & spray pressure
monitoring ensuring a clean wash

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
MIEL-PG8593

VALID UNTIL:
30 June 2024