

# Sliding Door Pharmaceutical Refrigerators

Includes **FREE** battery backed power failure and remote alarm box

Two sliding doors ensure maximum capacity with minimal footprint and forced air circulation for excellent temperature uniformity.

## Common Features:

- Temperature control range: 2°C to 14°C
- Stainless steel interior
- Alarms for high and low temperature, door open and power failure (with remote alarm box fitted)
- One access port
- Castors for mobility & locating feet when in working position
- Keyed door locks for security
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5-year warranty



**MPR-S163-PE, 158L**  
with 2 clear coated wire shelves

**Only \$2,995**

**Usually \$4,748**

Code: SANYMPR-S163-PE

- Cyclical (and forced) defrost
- LED interior lighting
- Quiet in operation: 35dB(A)



**MPR-514R-PE, 486L**  
with 5 polyester coated wire shelves & 5 roll out sliding racks

**Only \$7,155**

**Usually \$11,008**

Code: SANYMPR-514R-PE

- Fully automatic forced defrost
- Fluorescent interior lighting
- Quiet in operation: 42dB(A)

*Pricing valid for current stock only*

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: **VALID UNTIL:**  
SANY-Pharma **30 June 2024**