

Airstream® Gen 3 Horizontal Laminar Flow Cabinet 1.2m

One of the quietest (only 55.8dB(A)) and safest cabinet choices for sample protection requirements.

Only \$6,975

Normally \$7,671

Code: ESCOLHG-4AG-F8

Includes UV light, front cover and 2 power outlets


Applications:

Plant Tissue Culture	Electronics Inspection
Media Plate Preparation	Medical Device Assembly



- Fixed height polycarbonate sash and UV resistant glass side panels allow entry of "natural" light
- ULPA filter (10x filtration efficiency than of HEPA filter) delivers ISO Class 3 air cleanliness
- Sentinal™ Gold microprocessor controller with digital display showing cabinet status
- Stainless steel spill retaining worksurface for ease of cleaning
- Stable and self-compensating airflow
- Real time airflow monitoring with alarms
- Standby mode provides more energy savings and maintains an ISO Class 5 work zone
- ISOCIDE™ Powder Coat – painted surfaces have silver ion impregnated coat to inhibit microbial growth
- 3-year warranty

Optional adjustable height stand available for only **\$895** (usually \$1,515) just add code ESCOSTC-4AO



Energy efficient DC ECM Motor

(with 70% energy savings compared to AC motor)

price valid for current stock only

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: VALID UNTIL:
ESCO-Airstream 31 July 2024