

Guardian 5000 Hotplate Stirrer

Modernize your lab with the Guardian Series Hotplate Stirrers designed with outstanding safety & superior performance in mind.

Only \$990
Usually \$1,645
Code: OHAU30500521

The Guardian 5000 G51HS07C Hotplate Stirrer is an essential piece of equipment, which will shine in any heating/stirring application.

- Reliable heating (Ambient +5°C – 500°C)
- Digital control (Backlit LCD Display)
- Smart safety features (Hot top indicator light, audible alarm, SafetyHeat™ shut off)
- Exclusive SmartHousing - chemical resistant with an innovative, easy-to-clean design that channels spills away from internal components and the angled control panel
- Versatile (mixes up to 15L with speeds of 60 rpm – 1,600 rpm)
- Designed for everyday use (178 mm x 178 mm ceramic plate)
- 2-year warranty



SafetyHeat



SmartHousing



Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
OHAU-G51HS07C

VALID UNTIL:
30 Sep 2024