

MATOS® Aria Eco Fridge

Only \$3,995

Usually \$4,515

Code: ROLLMATOSARIA300R



Affordable, Compact & Reliable Cold Storage

The MATOS® Aria Eco fridge series is a new range of cost-effective European-built Refrigerators featuring the latest technology in controllers and compressors for accurate temperature stability.

With open door alarm, lock, and double glazed, toughened and insulated safety-glass doors.

Additional features include:

- Operating temp range: 1°C to +10°C
- Capacity: 269L Net / 290L Gross
- External Dimensions (mm): 597w x 645d x 1700h
- 6 x fitted shelves
- Controller: External microprocessor with LCD display
- 1-year warranty (extended 3-year warranty available)

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code:
ROLL-Aria300R

VALID UNTIL:
30 June 2024