

## Upright -86°C Freezer

### Proven Reliability

Large Capacity -86°C ultra-low temperature upright freezers with cascade refrigeration systems, available in both 728 and 845 litre capacities.

**Only \$14,900**

**Usually \$22,945**

Code: SANYMDF-U74V-PE



Storage capacity:  
728 litres = 528 x 2" cryoboxes  
Quiet in operation: 51dB(A)

**Only \$18,700**

**Usually \$28,791**

Code: SANYMDF-DU900V-PE



Storage capacity:  
845 litres = 616 x 2" cryoboxes  
Quiet in operation: 52dB(A)

**In Stock**

- Cascade refrigeration with Cool Safe compressors – proven reliable cooling technology with compressors running at lower temperatures.
- Two insulated gasketed inner doors with eagle-latches reduces heat loss.
- Temperature control range: -50°C to -86°C
- Hot gas bypass combined with fivefold door gasket assists in minimising frosting and ice build up around door seal
- Alarms for high and low temperature, door open, filter and power failure
- Three adjustable stainless steel shelves

- Three access ports
- VIP Plus insulation: provides typically 24% more storage capacity compared to conventionally insulated freezers with the same footprint
- Castors for mobility & locating feet when in working position
- Programmable start delay (in the event of a power failure to prevent overloading)
- Keyed door locks for security
- Dry remote alarm contacts
- Suitable for operation on 220/240V 50Hz 10 amp supply
- 5-year warranty

*Pricing valid for current stock only*

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.

Ref. Code: SANY-Upright-86  
VALID UNTIL: 31 July 2024