

Only \$6,925

Usually \$9,261

Code: GRANECOCOOL-100R

General Purpose 4.5L

Only \$7,995

Usually \$10,694

Code: GRANECOCOOL-150R

Advanced 5.5L



Innovative, eco-friendly, refrigerated and heating circulating baths offering significant running cost savings whilst delivering powerful cooling.

Consisting of two models, all products in the ecocool range are supplied assembled as ready to use kits, complete with accessory hosing, clips and connectors as standard.

- Temperature range -30°C to 150°C (model dependent)
- Active cooling through the whole temperature range
- True energy saving of up to 80% against standard compressor units
- Built in Britain to the highest specifications
- Thermostat and chiller work in harmony, neither will operate alone, eliminating any danger of overheating or freezing
- Single front switch for user convenience
- Modern, sleek, attractive design
- 3-year warranty
(with Grant renowned service and support, no registration required)

Applications

- Pharmaceutical - mini pilot plant reactors
- Education - rotary evaporator cooling, replacement of running tap water cooling, immersing small samples, photometry, chromatography systems
- Industrial - QC testing, sample preparation, general cooling, reaction chemistry, temperature control, semi-conductor manufacturing, rheometry
- Food - refractometry
- Life-science - electrophoresis cooling