

Miele

PG8583 Under Counter Washer/Disinfector

Only \$14,726

Was \$16,363

Code: MIEL62858306

accessories not included

Remove the chance of glassware breakage and injury during manual washing.

- Passive EcoDry drying function – activates as temperature descends to 70°C opening the door slightly enabling the load to be dried
- Two wash levels
- Spray arm and spray pressure monitoring ensuring a clean wash
- Stainless steel interior/exterior for durability and cleanliness
- Powerful 500L/min pump with integrated heating elements for safe washing cycles up to 93°C
- ADP model including pump for AD water and waterproofing kit
- Dispensing pumps for liquid detergent and liquid neutraliser
- 2-year Warranty



With many options and accessories to choose from contact us to assist in configuring your washing solution.



Ref. Code:
MIEL-PG8583

VALID UNTIL:
31 Oct 2024



0800 34 24 66

sales.nz@bio-strategy.com

www.bio-strategy.com

shop.bio-strategy.com

Prices and discounts may be subject to change and only available while stocks last. Not valid with any other special offer or discounts and cannot be applied in retrospect. All prices in local currency excluding GST and local delivery charges. Products advertised may be for research use only (RUO) or in vitro diagnostic use (IVD). Please refer to manufacturer's product insert for further information. For product use outside of that specified on the manufacturer's labelling and/or instructions for use, Bio-Strategy does not take responsibility.