

# NordicSafe® Class MSC (Low Noise)



- Equipped with energy-efficient ECM blower from ebm-papst Germany
- Low energy consumption (190 watts) and quiet

## ESCO NordicSafe Class II Microbiological Safety Cabinet: 1.2m with motorised sash

The Industry's  
Most Comfortable  
Energy Efficient Cabinet

Certified by PHE (Public Health England)  
for conformance to EN 12469

**Only \$14,694**

**Usually \$19,592**

Code: ESCONC2-4L8

Including **FREE**

UV light & 2 power sockets



Optional adjustable height stand for  
**\$950** (usually \$1,425)  
just add code ESCOSTC-4AO

### Operator protection, product protection and environment protection!

- Motorised sash ensures sash travels to working height automatically
- Ergonomic sloping front combined with UV resistant tempered glass side panels allow entry of more natural light
- ULPA downflow and exhaust filters with 99.999% efficiency, 10X better than HEPA filter (99.99%), improves safety and delivers ISO Class 3 air cleanliness
- Sentinel Gold microprocessor controller with digital display shows cabinet status
- Audible and visual alarms
- Airflow sensor monitors real-time airflow for safety and alerts the user if airflow is insufficient
- Night set back mode reduces power consumption by a further 60%
- 4-year warranty